

A GMP-compatible process for VSV-G+ and VSV-G-deficient LVs

The current 3rd generation lentivectors (LVs) bioproduction process developed at Alaya.bio is GMP-compatible for all types of vectors:


- Intracellular & surface expressed Genes of Interest (GOIs) with transgene sizes ranging from 2 to 8.9 kb
- With or without VSV-G envelope
- Integrative or non-integrative
- Mono- or bicistronic GOIs
- Fluorescent LVs (p17-GFP fusion)

✓ USP process in suspension cells can be performed in bioreactors




VSV-G pseudotyped lentiviral vector

- VSV-G envelope protein which confers broad tropism
- Produced by transfection of four plasmids
- Low transduction in non-activated/non-dividing-cells




VSV-G-deficient lentiviral vector

- Envelope protein: none
- Produced by transfection of three plasmids
- Deficient for transduction




Targeting nanoparticles (TNPs) platform for *in situ* gene delivery




QCs of 24 R&D batches at the 2 L scale


Viral particle content of LV batches



Functionality of LV batches



Process-related impurities of LV batches



At the end of the entire process, 20 to 30 mL of sterile purified and filtered drug product are obtained (one point = one batch).


- LV product quality is monitored based on:
- ✓ Viral particle concentrations (physical, p24-based and RNA-based titers)
 - ✓ Infectious titers determined by qPCR on Jurkat cells (3 days post-transduction)
 - ✓ Process-related impurities contents (Host Cell Proteins and residual DNA)

Downscaling for the evaluation of new constructs

USP process performance from 18 mL to 2 L

USP process parameters routinely operated for the manufacturing of LV DPs for our own developments were tested and transferred to the 30 mL- then the 18 mL- scale to develop a screening tool for new LV variants based on their manufacturability at the 2 L scale and their functionality.

- Physical titers were measured by interferometry, anti-p24 ELISA and RT-qPCR
- Functionality was assessed for VSV-G+ vectors by transduction of Jurkat cells




- Data shown here highlight the robustness of the USP process with VSV-G+ and VSV-G-deficient CD19-CAR & GFP-encoding, non-purified and clarified bulk LVs (one bulk = one point)
- Variants of the CAR construct (such as CD19-CAR-P2A-Reporter gene, CD19-CAR & p17-GFP fusion) were engineered and the feasibility of their bioproduction was evaluated. Consistent titers from the 18 mL- up to 2 L-scale validated the power of the downscaling screening approach
- ✓ The predictability of yields at the 18 mL-scale now allows us to evaluate, at lower cost and in a shorter timeframe, new genetic constructs provided by potential partners

Downscaling to improve yields and scalability of DSP process units

- Challenge: Low recovery yields and high batch-to-batch variability during clarification and concentration-diafiltration steps
- Identified potential root cause: oversized surface area
- Small Scale Filtration Systems Models to lower material consumption, increase throughput of tested conditions and rapidly identify most efficient scale


Case study #1: Clarification optimization by depth filters sizing model



The downscaling of the clarification step has been validated by 3 independent studies based on process parameters (pressure, turbidity) and LV integrity (RNA-and p24-based titers) and functionality (infectious titer).

- ✓ P_{max}™, V_{max}™ and T_{max}™ sizing methods were tested and identified 2-fold smaller filtration areas with comparable performances
- ✓ Best small-scale parameters have been successfully validated with our current 2 L-scale bioprocess, thereby confirming the scalability of such optimization
- ✓ An additional bioburden reduction filter has been implemented to tightly control the quality of filtrates
- ✓ A turbidity extrapolated-model was established to anticipate any rise of Bulk impurities and size the most suitable depth filter train compatible with scale-up

Case study #2: TFF optimization by cut-off comparison with Ambr® Crossflow



A TFF scale-down model was conducted with our current 2 L-scale bioprocess parameters (C/D/C mode, pressures, diavolumes and diafiltration buffer) and analytical read-outs. Two loading volumes were tested and revealed different performances between 100 and 300 kDa 50 cm² membranes:

- ✓ Oversizing of current cassette area (300 kDa - 0.14 m²) was detrimental for particle recovery due to LV entrapment in pore mesh (at 10 L/m² - our reference parameter) despite a reduced process time
- ✓ Small scale model revealed that a 2-fold reduction of surface area may translate into a 2-fold gain in recovery yields
- ✓ Two options now available: 10 L/m² & 100 kDa membrane with wash step = recovery x4 vs 40 L/m² & 300 kDa membrane no wash = recovery x3
- After Filtration area optimization, loading pressure is under evaluation

Cross-functional analytical tools

Complementary Quality Control methods are used for in-process control of the USP/DSP process performance and the evaluation of the quality of the vector as final drug product or payload for TNPs.


Quality attributes of LV payloads

- Viral p24 capsid protein quantification (ELISA)
- Proviral viral RNA copies quantification (RT-qPCR)
- Particle concentration, size and aggregation (NTA, interferometry)
- Infectious titers (qPCR post-transduction) & GOI functionality (cytometry)
- Zeta potential and viral particle diameter (DLS)

Process-related impurities

- Host Cell Proteins quantification (ELISA)
- Residual dsDNA quantification (Quant-IT™ PicoGreen™)

Alaya.bio platform highlights



- ✓ A universal and GMP-compatible bioprocess developed for the production of high-quality LV batches for ex vivo or *in situ* applications
- ✓ Cost- and time-effective evaluation of the manufacturability of new constructions for internal or third-parties is achieved with our small-scale production platform
- ✓ Our bioproduction platform is amenable to non-integrative LV alternatives when looking for transient expression of therapeutically relevant GOIs

Got questions or just curious?
We'd love to hear from you - hello@alaya.bio