

AudioNova V

Product Overview 3.0





AudioNova V offers our best sound performance yet

What's new with AudioNova V?

SoundSuite OS™ now includes a new listening environment - conversation in loud noise - for a total of 8 different listening situations.

ActiveFocus improves SNR for one-to-one conversations in noisy environments

The new 360 conversation in car manual program focuses on speech automatically for passengers, regardless of where they are seated.

SoundSuite OS

AudioNova V features SoundSuite OS™, the highly advanced signal processing system that delivers an adaptable and personalized solution.

MFA Connectivity

Patients can make hands-free phone calls and stream music, podcasts, and other media directly into their hearing instruments from any Bluetooth® enabled device with Made-For-All capabilities. AudioNova V can even connect to two different devices at once! And with tap control, they can control media and calls with a quick tap to the ear.

Hearing Remote App


Your patients can personalize their listening experience by improving focus, boosting speech, and reducing background noise using the Hearing Remote App*. The app also allows them to select preset manual programs that suit their specific lifestyles. And for changes that require your expertise, remote adjust allows you to deliver professional fine-tuning remotely. This gives flexibility and convenience to all your patients by facilitating easy adjustments without an in-office visit.

AudioNova:upgrade

Your patients can upgrade their hearing solution, thanks to AudioNova:upgrade. Patients pay the difference between their current and upgraded technology level, without the need to purchase new hearing aids. AudioNova:upgrade can help you overcome initial patient cost concerns, by starting them with a technology level that is comfortable and affordable to them with the ability to upgrade as their hearing needs change.

AudioNova V Family

AudioNova V has new form factors, that focus on aesthetics, comfort and ease-of-use.



Behind-the-ear (BTE)

Receiver-in-canal (RIC)



AudioNova V-UP

Our Ultra Power BTE with a traditional battery and direct connectivity



AudioNova V-PR

Our rechargeable BTE with direct connectivity



AudioNova V-R

Our smallest rechargeable RIC with direct connectivity



AudioNova V-RT

Our rechargeable RIC with telecoil and direct connectivity



AudioNova V-312

Our smallest RIC with a traditional battery and direct connectivity

Accessories

RIC chargers

The RIC chargers feature magnetic retention, allowing your patients to remove their devices with ease.



V charger for AudioNova V-R/RT

BTE chargers

Charger can be combined with optional Power Pack.



Easy Line BTE charger for AudioNova V-PR



Easy Line Combi BTE2 charger for AudioNova V-PR



TV Connector

Your patient's favorite media streamed wirelessly in stereo sound with TV Connector.




PartnerMic

Is an easy-to-use lapel microphone ideal for one-on-one conversations.



Remote Control

Your patient can discreetly adjust volume and switch programs with the Remote Control.



Hearing Remote App

The Hearing Remote App works on iOS and Android smartphones. It puts personalization at your patient's fingertips.

At a glance




Fitting guides


SoundSuite OS™

		Level 9 Premium	Level 7 Advanced	Level 5 Standard	Level 3 Essential	Level 1** Lower Essential
Environmental classification	Total listening environments	8	6	4	2	AutoMic
	Noise	•				
	Conversation in loud noise	•				
	Conversation in noise	•	•			
	Music	•	•			
	Quiet	•	•	•		
	Conversation in a small group	•	•	•		
	Conversation in a crowd	•	•	•	•	
	Conversation in quiet	•	•	•	•	
	Total streaming environments	2	2	2	2	2
Media speech	•	•	•	•	•	
Media music	•	•	•	•	•	
Sound optimization and performance	Auto Target Pro	•				
	Auto Target		•			
	ActiveFocus	•	•			
	Dynamic noise reduction	•	•			
	Soft speech intensifier	•	•	•		
	Sound Mapping	•	•	•		
	Pinna Effect	•	•	•	•	
	Speech intensifier	•	•	•	•	
	Noise reduction	•	•	•	•	
	Hardware features	Tap control*	•	•	•	•
Made-For-All direct connectivity		•	•	•	•	
Rechargeable		•	•	•	•	
Telecoil*		•	•	•	•	
Fine-tuning channels		20	20	16	12	


RIC




S receiver



M receiver



P receiver



UP receiver

- Open dome/cap dome
- • • Vented dome
- Power dome or sleeve mold

Available in all technology levels

Sound stabilization

- Pulse Protector 2
- Wind control
- Feedback manager
- Natural Sound

Ease and convenience

- Hearing Remote App
- TV Connector
- Remote Control
- PartnerMic
- RogerDirect
- Wireless synchronization
- Binaural Phone[‡]

Microphone options

- Adaptive directional
- Fixed directional
- Fixed wide directional
- Omni directional

Fitting

- First fit approach
- Automatic Adaptation Manager
- Frequency compression 2
- Tinnitus masker
- Manual programs
- IntelliVent

Experience Innovations

- AudioNova:trial
- AudioNova:upgrade

Digital Solutions

- Tutor
- Remote adjust
- Capture All
- Data logging

Personalization

- Equalizer
- Clarity and comfort buttons[‡]
- Optional app programs[‡]

*Compatibility may be limited by product style.

**Not all products in all markets available at technology level 1.

‡ Not available for technology level 1.